

## ZETASIZER CELLS AND CUVETTES: DESCRIPTIONS AND USES

A wide selection of cells is available for use with the [Zetasizer range](#) of instruments. This technical note lists all of the cells currently available from Malvern Panalytical. Each has qualities and advantages that suit particular applications. The following tables describe the available cells and show which instruments from the Zetasizer range that the cells can be used with.



Table 1 gives a description of the cells.



Table 2 shows which of the cells are compatible with size measurements.



Table 3 shows which of the cells are compatible with zeta potential measurements.

Table 4 shows which of the cells are compatible with molecular weight measurements.

Table 1: Size, molecular weight and zeta potential cells available with the Zetasizer range.



Part number	Description	Picture
DTS0012	Disposable polystyrene cuvette, pack of 100 with 100 caps, minimum volume 1mL, suitable for all systems that measure size. Compatible with: Zetasizer 1000/3000, HPPS, Zetasizer Nano series	
ZEN0118	Disposable low volume cuvette for size measurement, pack of 100 with 100 caps, minimum volume 50µL, suitable for all systems that measure size at 90 degree scattering angles only. Can be used with aqueous dispersants. Compatible with: Zetasizer Nano series and Zetasizer µV	


Part number	Description	Picture
ZEN0040	<p>Disposable, solvent resistant micro cuvette, for size measurement at a 173 degree scattering angle, pack of 100 with 100 caps, minimum volume 40<math>\mu</math>L, suitable for Zetasizer Nano S and ZS. Not suitable for measurements at a 90 degree scattering angle. (Resistant to Acetone, Benzaldehyde, Butanone, Dioxane, DMF, Ethyl acetate, Isopropanol, various acids and bases)</p>	
PCS8501	<p>Square glass cuvette for size and molecular weight measurements. For the Zetasizer 1000/3000/3000HS, HPPS, Nano series and Zetasizer <math>\mu</math>V. With circular aperture and stopper. Not suitable for use with the 'dip' cells for electrophoresis measurement</p>	

Part number	Description	Picture
PCS1115	<p>Square glass cuvette with cap for 'dip cell' electrophoresis measurement and size and molecular weight measurements. Compatible with the Nano series universal dip cell (ZEN1002) and the Zetasizer 2000 and 3000 Series aqueous and non-aqueous 'dip' cells (DTS5001/5002). With square aperture and cap</p>	
ZEN0023	<p>Quartz flow cell for size, intensity and molecular weight measurement. Compatible with MPT-2 autotitrator and Zetasizer Nano S, ZS, S90. (S90 for Size and Intensity only). Complete with connecting tubing.</p>	

Part number	Description	Picture
ZEN0116	<p>Quartz flow cell to enable the use of the Zetasizer Nano as a detector when connected to an SEC system. Connections compatible with standard SEC tubing 1/16inch o.d. Includes, ZEN0023 quartz flow cell, connections for cell, 1.5m PEEK tubing, 1/16inch o.d. 0.1mm i.d.</p>	 <p>A photograph of the ZEN0116 quartz flow cell. It is a black, rectangular device with two clear quartz tubes protruding from the top. The front face features a label with 'QS' in a blue box, '3.00 mm' below it, and a white downward-pointing arrow. A small white rectangular window is visible below the arrow.</p>
ZMV1008	<p>Quartz flow cell to enable the use of the Zetasizer <math>\mu</math>V as a detector when connected to an SEC system. 8<math>\mu</math>L volume. Connections compatible with standard SEC tubing 1/16inch o.d. Includes, quartz flow cell, connections for cell, 1.5m PEEK tubing, 1/16inch o.d. 0.1mm i.d.</p>	 <p>A photograph of the ZMV1008 quartz flow cell. It is a black, rectangular device with two clear quartz tubes protruding from the top. The front face features a label with 'QS' in a blue box and a white downward-pointing arrow. A small white rectangular window is visible below the arrow.</p>

Part number	Description	Picture
ZEN2112	<p>Low-volume quartz batch cuvette for size and molecular weight measurements for use with the Zetasizer Nano ZSP, ZS and S. Minimum sample volume is 12µL.</p>	
ZEN1002	<p>Universal Dip Cell for the zeta potential measurement of low mobility samples such as those in organic solvents. Used in conjunction with DTS0012 and PCS1115. Compatible with the Zetasizer Nano ZSP, ZS, Z and ZS90.</p>	

Part number	Description	Picture
ZMV1002	Low-volume quartz batch cuvette for size and molecular weight measurements for use with the Zetasizer $\mu$ V. Minimum sample volume 2 $\mu$ L	
ZMV1012	Low-volume quartz batch cuvette for size and molecular weight measurements for use with the Zetasizer $\mu$ V. Minimum sample volume 12 $\mu$ L.	

Part number	Description	Picture
ZEN1020	<p>Surface zeta potential cell for the measurement of the zeta potential of planar surfaces using the Zetasizer Nano ZSP, ZS, Z and ZS90. Used in conjunction with DTS0012 and PCS1115. Samples should be 4-7 mm x 4 mm x <math>\leq 1.5</math> mm thick.</p>	 <p>The image shows a ZEN1020 surface zeta potential cell. It is a vertical, cylindrical device with a black top cap and a clear quartz body. A silver-colored metal band is visible around the middle section. The top cap has a small window showing a blue and green display. The device is designed for measuring the zeta potential of planar surfaces.</p>
ZEN1010	<p>High concentration cell for the measurement of size in the Zetasizer Nano ZSP and ZS, and the measurement of zeta potential in the Zetasizer Nano ZSP, ZS and Z. The quartz cell has a path length of 1.5mm and this cell is compatible with the MPT-2 autotitrator.</p>	 <p>The image shows a ZEN1010 high concentration cell. It is a vertical, rectangular device with a black top cap and a clear quartz body. The device is designed for measuring the size and zeta potential of high concentration samples. It is compatible with the MPT-2 autotitrator.</p>




Part number	Description	Picture
DTS1070	<p>Completely disposable cuvettes for the measurement of zeta potential with the Zetasizer Nano ZSP, ZS, Z and ZS90. It is also possible to measure size in these cuvettes with the Zetasizer Nano ZSP, ZS and S. These cuvettes are compatible with the diffusion barrier technique of sample introduction and the MPT-2 autotitrator. Pack of 10 with 25 stoppers.</p> <p>Requires version 7.02 or later of the Zetasizer software.</p>	

Table 2: Cells capable of size measurements with the different Zetasizers.

Cuvettes for size measurement	Zetasizer 1000/3000	HPPS	Zetasizer Nano S90/ZS90	Zetasizer Nano S/ZS/ZSP	Zetasizer $\mu$ V
DTS0012	✓	✓	✓	✓	✓
ZEN0118	✓	✗	✓	✗	✓
ZEN0040	✗	✗	✗	✓	✗
PCS8501	✓	✓	✓	✓	✓
PCS1115	✓	✓	✓	✓	✓
ZEN0023	✗	✗	✓	✓	✗
ZEN0116	✗	✗	✓	✓	✗
ZMV1008	✗	✗	✓	✗	✓
ZEN2112	✗	✗	✗	✓	✗
ZMV1002	✗	✗	✗	✗	✓
ZMV1012	✗	✗	✗	✗	✓

Table 3: Cells capable of zeta potential measurements with the different Zetasizers.

Cuvettes for zeta potential measurement	Zetasizer Nano Z	Zetasizer Nano ZS90	Zetasizer Nano ZSP/ZS
DTS1070	✓	✓	✓
ZEN1020	✓	✓	✓
ZEN1010	✓	✗	✓
ZEN1002	✓	✓	✓

Table 4: Cells capable of molecular weight measurements with the different Zetasizers.

Cuvettes for molecular weight measurement	Zetasizer Nano ZS90	Zetasizer Nano ZSP/ZS	Zetasizer $\mu$ V
PCS8501	✓	✓	✓
PCS1115	✓	✓	✓
ZEN2112	✓	✓	✓
ZMV1002	✗	✗	✓
ZMV1012	✗	✗	✓
ZMV1008	✗	✗	✓
ZEN0116	✓	✓	✓

## MALVERN PANALYTICAL

Groewood Road, Malvern  
 Worcestershire, WR14 1XZ  
 United Kingdom  
 Tel. +44 1684 892456  
 Fax. +44 1684 892789

Lelyweg 1,  
 7602 EA Almelo,  
 The Netherlands  
 Tel. +31 546 534 444  
 Fax. +31 546 534 598

info@malvernpanalytical.com  
[www.malvernpanalytical.com](http://www.malvernpanalytical.com)

Disclaimer: Although diligent care has been used to ensure that the information in this material is accurate, nothing herein can be construed to imply any representation or warranty, as to the accuracy, correctness or completeness of this information and we shall not be liable for errors contained herein or for damages in connection with the use of this material. Malvern Panalytical reserves the right to change the content in this material at any time without notice.  
 Copyright: © 2025 Malvern Panalytical. This publication or any portion thereof may not be copied or transmitted without our express written permission. [tn120917zetasizercells.cuvettes-EN](https://www.malvernpanalytical.com)